

Start-up system

- Switch on the power socket located on the floor
- Switch devices on in the order of:
 - Mercury lamp control-unit switch on and ignite
 - Halogen lamp
 - Incubator control-box
 - Microscope stand located at the right, behind the serial plug
 - Stage controller located at the back of the device
 - Shutter controller
 - Camera via power supply located on the floor
- Wait approximately 5 seconds before you start the computer
 - Log in with user: nikon01 password: nikon01
- Start imaging software Metamorph
 - Select "User" account

Using Tokai Hit Incubator + CO₂

- Fill reservoir with demineralized water up to level of the CO₂ tube
- Place lens heater around desired objective (adjust the Velcro to fit around lens)
 - For Oil lenses only
- Settings for the heating unit are:

	Dry objective	Oil objective
○ Top heater	39.0 °C	41.5 °C
○ Bath heater	36.5°C	38.5 ° C
○ Stage heater	36.5°C	38.5 ° C
○ Lens heater	-	37.0 ° C
- Open the valves of the CO₂ gas cylinder
 - Turn the valve on top of the cylinder ¼ open
 - Turn the tap at the pressure meters open
 - Pressure at the left meter should be 1-2 bar
- Enable on CO₂ mix unit of the incubator control-unit
- Percentage is set to 5%, mixing gasses takes about 5-10 minutes

- Elements of the incubator can be switched on/off individually
 - Use this is when lens heating is not needed
- Allow the system to warm up for ~30 minutes

Focus your cells

- Choose your objective
 - 100x , 60x , 40x **oil** or 20x **air**
- Put the objective in the lowest position
 - By using the focus knob, or button Escape (right side lowest button)
 - Lowest position is ~500 um
- Put a little drop emersion oil on the lens
 - Only if it's an oil emersion objective
- Centre you coverslip/well/cells above the objective with the joystick
 - Button on top of joystick will change speed, S (slow) or F (fast)
- Focus the objective till the oil is in contact with the coverslip
- In Metamorph for normal halogen light choose:
 - for 20x: Taskbar button: DIC Binocular
 - use filter "N2" in microscope head
 - for 100/60/40x: Taskbar button: DIA Binocular
 - use filter "A" in microscope head
- Adjust precision of focus wheel
 - course, fine or extra fine Buttons are located left and right on the stand
- The display of the microscope shows the height of the objective
- Buttons L100 and EYE at the front of the microscope stand display the direction of the light path (left -> camera, up -> eyepiece)

Hardware autofocus: Perfect Focus System (PFS)

- When imaging with just a coverslip, hardware autofocus is possible
 - Orange LED will light up if coverslip is within reach
 - Enable PFS with button "Focus" located at the front of the microscope stand
 - Focussing is done with external focus device
 - Button at the front changes precision of focus wheel
 - Focusing with PFS should be viewed through the camera, it's to precise to see with your eyes

Controlling intensity mercury lamp

- Fluorescence can be detected with use of the Mercury lamp
 - Choose illumination settings according to your fluorophore or dye
 - By pressing a button in the taskbar
- Intensity of the light can be controlled with neutral density (ND) filters.
- These are located between the actual lamp (right) and the microscope stand.
 - Press inwards to reduce intensity by 4 or 8 times (ND4, ND8)
 - Most right filter will block all light

Metamorph - Taskbar


- Metamorph functions from a taskbar and an Acquire window
- When the taskbar is not visible
 - Load the taskbar via: Journals>Taskbar>Load Taskbar
 - Location is C:\MM\TASKBARS\
- You can reset the illumination settings with the first button: "Restore OIC settings"
- Select preferred illumination with buttons
 - Left column to eyepiece
 - Right column to camera with live mode

Restore illumination settings	Main Taskbar (OIC)	Open acquire window
Open stream window		Open MDA configuration
Select illumination setting to <u>eyepiece</u>		Select illumination setting to <u>camera</u> and start Live mode

Metamorph - Acquire window


Sidebar left:

- Exposure time
 - Minimum interval time at full chip is 90ms, quad chip is 50ms
 - Exposure time can be lower than interval time
- Full Chip, Quad Chip, Use Active Region
 - Use full chip, middle quadrant, or region selected by ROI
- Live bin & bin
 - Combine pixels to increase readout speed
 - Use only for live viewing or real recording
 - Normal state = 1
- Start Live/Stop Live
 - Start/stop live mode




Acquire images

- To acquire live mode press "Acquire" button
- Separate images can be recorded as separate images or in stack



Tab: Display:

- Auto scale option
 - Enable Auto scale or change custom auto scale
 - Full bit-range is saved to file, this is only for visualization in the software
- Visualize saturated pixels




Tab Acquire:

- Be sure to set the Illumination choice to "Current shutter"


Tab Special:

- Adjust digital gain
 - Calibrated gain: 1, 2, 3
 - Gain 3 (4x) is standard option
- Use Visual Gain to control gain between 0 and 580
- Use "Frames To Avg" option to average images / or use longer exposure time



Live image window

- Magnifying glass
- Display image in other LUT
 - Grayscale, Colour by wavelength
- Histogram of intensities
 - Manual scaling can be done with the orange arrows
- Auto scaling options
 - Best fit range, 8/12/14/16 bit
 - Auto scale on/off
 - Scale within region



Streaming imaging

Open the Stream Acquisition window with the button in the


- Start recording with Acquire button
- Check if chosen configuration is ok

Tab Acquire:

- Acquisition modes:
 - RAM save later to hard disk
 - HD choose name and location to save
Save as on file or separate TIF files
- Enter the desired number of frames
- Select "Initial illumination": *Current Shutter*

Tab Camera Parameters:


- Recording options:
 - Average frames
 - Number of frames to skip, during acquisition
 - Don't show recording (faster when exposure time < minimal interval time)
- Camera readout default options should be:
 - Camera state: HALT
 - Shutter mode: OPEN PRE SEQUENCE
 - Clear mode: CLEAR PRE SEQUENCE



Multi-Dimensional Acquisition (MDA)

Combine multiple tasks: z-stack, timelapse, streaming, multiple positions, multiple wavelengths, use of journals.

- Open the MDA window using Taskbar
- Select acquisition mode(s) in *Main*
- Continue with the tabs shown left:
 - Saving: select filename and folder
 - Select illumination settings
 - ... options for acquisition modes
- When setup is done press acquire
- Save or load settings of experiment setup, with Save/Load State buttons
- Live mode can be started from MDA window



MDA: Time lapse imaging

- Activate the time-lapse option in MDA's main window
- Set the desired interval and time
- Choose the desired illumination settings in the wavelength tab
 - Or multiple when the option multiple wavelengths are selected
- Start recording by clicking on Acquire

MDA: Z-stack

- Activate Z-stack option in MDA's main window
 - Activate use dual Z motor
 - Set all options to: *TiZ*
- Select option: Center around current
 - Height of objective is shown (same as display microscope)
- Give range, step size (in um), number of steps
- Choose the desired illumination settings in the wavelength tab
- Start recording by clicking on Acquire
- Don't use PFS during Z-stack
 - Activate "Run Journal" in MDA's main window
 - Choose PSF off before Z-stack and PFS on after Z-stack

Scan Slide (a.k.a. Tile scan)

Use this to make images of a larger area and stitch them together

- Access via Apps > Scan Slide
- Main:**
- Load Settings for re-use purposes
 - Load calibration for used magnification

Acquisition:

- Select and configure desired wavelength(s)


Slide area:

- Navigate to upper left region and click on "set to current"
- Navigate to lower right region and repeat

Data review:

All images are saved, last image in sequence is stitched image.

- Stitch file in Data review tab
 - Choose overlap and channels



Shut down system

- Save data to OIC network storage (O:\ drive)
 - Connect to network storage with your own PC via the address: <\\oic-station\oic>
user: guest, password: guest
 - When the hard disk of the microscopes pc gets full, data will be deleted
 - You are responsible for your own data
- For live cell imaging:
 - Remove Tokai Hit incubator from the stage and place it on the table
 - Close the Tokai Hit incubator with a plastic well
- **Clean objective carefully**
 - Remove oil with lens paper:
 - Hold lens paper tight between your fingers and wipe over objective
 - Clean lens with 2-propanol on your lens paper
- Cover stage with lens paper cleaning box
- **If you are the last user of the day**
- Shut down all devices
 - Shut down Mercury lamp at last
 - Camera power is switched at the adapter box at the floor
- Switch off the power socket located at the floor

OIC contact information:

- Martijn de Gruiter Be-341a tel. 42205
- Gert-Jan Kremers Be-346 tel. 43578