


### Start-up system

- Switch on the power socket (#1)
- Switch on Mercury lamp control unit      switch on and ignite (#2)
  - Wait ~10 seconds for the device to load
  - Hold ignite button until orange LED is active
- **Start the PC before other hardware is active**
  - User account Nikon1078, no password required
- Switch on the power socket (#3)
- Switch devices on in the order of:
  - Halogen lamp
  - Incubator control-box
  - Microscope stand                      located at the right, behind the serial plug
  - Stage controller                        located at the back of the device
  - Shutter controller
- Start imaging software Metamorph
  - Select "User" account
- Upon error; drag error message from behind Metamorph logo
  - Message: *ASI MS2000 not detected on COM1*
  - Choose: *Retry*



### Shut down system


- For live cell imaging:
  - Remove Tokai Hit incubator from the stage and place it on the table
  - Close the Tokai Hit incubator with a plastic well
- **Clean objective carefully**
  - Remove oil with lens paper:
  - Hold lens paper tight between your fingers and wipe over objective
  - Wet your lens paper with 2-propanol and wipe 3-4 times more
- Cover stage with lens paper cleaning box
- **If you are the last user of the day**
- Shut down all devices
  - Shut down Mercury lamp at last
  - Camera power is switched at the adapter box at the floor
- Switch off the power socket located at the floor

### Controlling intensity mercury lamp

- Fluorescence can be detected with use of the Mercury lamp
  - Choose illumination settings according to your sample
- Intensity of the light can be controlled with neutral density (ND) filters.
  - Located between the actual lamp (right) and the microscope stand.
  - Press inwards to reduce intensity by 4 or 8 times (ND4, ND8)

### Taskbar


- When the taskbar is not visible
  - Load the taskbar via: Journals>Taskbar>Load Taskbar
  - Location C:\MM\TASKBARS\

Restore illumination settings		Open acquire window
Open stream window		Open MDA configuration
		Move stage with cursor
Select illumination setting to <u>eyepiece</u>		Select illumination setting to <u>camera</u> and start Live mode

- If needed: reset the illumination settings with: "Restore OIC settings"
- Open the Acquire window
- Configure Stream acquisition or start MDA configuration
- Move stage by clicking inside live image window (distance from centre equals movement)
- Select illumination settings for Eyepiece (left) or camera (right)

## Acquire window

- Open via Taskbar button “Acquire”
- Acquire an image
- Set exposure time
- Camera gain settings
- Show/Stop live view



## Multi-Dimensional Acquisition (MDA)

Combine multiple tasks: z-stack, timelapse, streaming, multiple positions, multiple wavelengths, use of journals.

Select acquisition mode(s) in *Main*

- Continue with the tabs shown left:
  - Saving: select filename and folder
  - Select illumination settings
  - ... options for acquisition modes
- When setup is done press acquire
- Save or load settings of experiment setup, with Save/Load State buttons
- Live mode can be started from MDA window

